# RA-915M

#### Zeeman mercury analyzer



# Mercury determination in natural gas according to ASTM D5954-98(2014)e1 and ISO 6978-2:2003

#### **INTRODUCTION**

ASTM D5954-98(2014)e1 and ISO 6978-2:2003 standards specify a method for determination of the mercury mass concentration down to 0.001 µg/m<sup>3</sup> in pipeline quality natural gas. These standard methods are based on collection of mercury by amalgamation on a gold/platinum alloy with the follow-up desorption of atomic mercury and its determination by means of AAS or AFS. Lumex Instruments has developed a practical guidance for implementing ASTM D5954-98(2014)e1 and ISO 6978-2:2003 methods using RA-915M mercury analyzer.

#### **MEASUREMENT METHOD**

The method includes two consecutive stages:

- Gas sampling: A known volume of natural gas passes through a sampler with installed sampling tubes and a flow meter; mercury is accumulated on the sorbent due to amalgamation.
- Analysis: The sampling tube is heated up to 700–800°C, the collected mercury is released and its mass is measured using RA-915M mercury analyzer.

# **MEASUREMENT RANGE**

The measurement range of the mercury mass concentration in natural gas is  $0.001-100 \ \mu g \ /m^3$ .

# **EQUIPMENT AND REAGENTS**

The following equipment and materials are used for analysis:

- Gas sampler for sampling tubes according to ASTM D5954-98(2014)e1 and ISO 6978-2:2003;
- Mercury analyzer RA-915M with PYRO-915+ and RP-92 attachments;
- PC with Windows® and RAPID software;
- Kit for mercury determination in natural gas (includes sampling tubes, CRM of mercury ions, and activated carbon with mercury content ≤2 ppb).

# **EXAMPLES OF ANALYSIS**





Reusable sampling tubes for mercury sampling from natural gas.



Procedure of the sorption trap analysis.

Results of analyses of the sampling tubes spiked with 0.5 ng of mercury: **1** – 0.514 ng; **2** – 0.515 ng; **3** – 0.506 ng. One measurement takes 2 minutes.



www.lumexinstruments.com | sales@lumexinstruments.com