



# DETERMINATION OF **INORGANIC ANIONS** IN AQUEOUS MATRICES ACCORDING TO EPA 6500 AND ASTM D6508-15 TEST METHODS

**EPA 6500 ASTM D6508-15** 

#### INTRODUCTION

These test methods are applicable for determination of the dissolved inorganic anions: fluoride, bromide, chloride, nitrite, nitrate, *ortho*-phosphate, and sulfate in aqueous matrices using capillary ion electrophoresis with indirect UV detection. The methods can be applied to **drinking water**, wastewater, and ground water.

## **MEASUREMENT METHOD**

The measurement method is based on capillary zone electrophoresis with indirect UV detection at the wavelength of 254 nm.

#### **MEASUREMENT RANGE**

Anions	Measurement range, mg/L
Chloride, bromide, nitrite, sulfate, nitrate, ortho-phosphate	0.1–50
Fluoride	0.1–25

#### **EQUIPMENT AND REAGENTS**

The CAPEL capillary electrophoresis system is used in measurements. Data acquisition, collection, processing and output are performed using a personal computer running under WINDOWS® XP/7/8/10 operating system with installed dedicated software package ELFORUN. All reagents must be of analytical grade or better.

Lumex Instruments set, order No 0300002708.

# **EXAMPLES OF REAL ANALYSES**


**BGE:** chromate, with TTAOH, CHES, and calcium gluconate


Voltage: - 15 kV

Temperature: 25 °C

Capillary:  $L_{eff}$ /  $L_{tot}$  50/60 cm, ID 75  $\mu m$  Detection: 254 nm, indirect


Injection: 300 mbar x sec





### **Measurement results:**

- 1 chloride (7.78 mg/L)
- 2 sulfate (26.3 mg/L)
- 3 nitrate (1.66 mg/L)
- 4 carbonate



The contents on this paper are subject to change without notice.

To get more specific information, please contact the representative by sales@lumexinstruments.com