

# RA-915 series

Zeeman AA spectrometer



## Total mercury determination in water by CVAAS

### INTRODUCTION

Mercury determination in surface, drinking and waste waters is one of the most popular analyses used for environment pollution and sanitary control.

### MEASUREMENT METHOD

Samples of surface, drinking and waste waters must be pre-treated using the appropriate standard digestion procedures (for instance, by potassium permanganate or bromide-bromate mixture).

The method is based on the reduction of Hg(II) to the atomic state using a tin(II) chloride solution and the follow-up transporting of mercury atoms into the analytical cell of the analyzer by air flow (the "Cold Vapor" technique). The mercury concentration is then measured by the Mercury Analyzer RA-915M combined with RP-92 or URP attachment. Mercury mass concentration is calculated by the analytical signal integration and pre-established calibration (area of the peak vs mercury mass).

Sample volume for injection is from 1 to 20 mL. Detection limit (DL) is 0.0005 ng of mercury, which corresponds to the mercury mass concentration of 0.000025 µg/L (0.025 ng/L). The measurement range is 0 to 5 µg/L for the multi-path cell and can be extended up to 2000 µg/L using the short single-path cell.



RA-915M + RP-92




Fig. 1. Peaks corresponding to the injection of 0.005 ng of Hg.

### EQUIPMENT AND REAGENTS

The following equipment and materials are used for analysis:

- Mercury analyzer RA-915M with RP-92 or URP attachment;
- PC with Windows® XP/7/8/10 and RAPID software;
- Lumex Instruments kit, **Order No 0300003071.**

### COMPATIBLE METHODS LIST

- EPA Method 245.1
- ISO 12846:2012
- ASTM D3223-17
- AOAC Official Method 977.22

Directives & standards for drinking water	Limits, µg/L
WHO Guidelines for drinking water quality (2011)	6
Drinking Water Directive 98/83/EC	1
US EPA National Secondary Drinking Water Regulations	2
TR EAEU 044/2017 Technical Regulation on Packaged Water	0.2 / 0.5 / 1
GB 5749-2006 Standards for drinking water quality	1




Fig. 2. Examples of analyses of the surface and waste water using EPA Method 245.1.

	Water samples	Measured value, µg/l
1	OK L4-14	0.51±0.10
2	OK 13-10	1.9±0.5
3	OK C3-13	5.2±0.7



Find your local Lumex Instruments distributor:  
[www.lumexinstruments.com](http://www.lumexinstruments.com)  
[sales@lumexinstruments.com](mailto:sales@lumexinstruments.com)