



## Determination of **lactose** in lactose-free milk processing products

### INTRODUCTION

Lumex Instruments proposed the method for determination of lactose content in **lactose-free milk processing products** by capillary electrophoresis. For the determination of lactose, fructose, glucose, and saccharose content in other types of products (including determination of lactose in products with lactose content of less than 1 %), it is recommended to use Lumex Instruments set **No. 0300002753**.



### MEASUREMENT METHOD

The method is based on extracting the component from samples with distilled water, obtaining derivatives, further separation, identification and determination of the mass concentration of lactose by capillary electrophoresis. The component is detected at a wavelength of 305 nm.

### MEASUREMENT RANGE

The measurement ranges for the mass concentration of lactose is **0.004 % – 0.7 % (0.04–7.0 g/L)**. Specially selected conditions for electrophoretic separation eliminate the interfering effects of glucose and galactose, which are present in significant quantities in lactose-free products.

### EQUIPMENT AND REAGENTS

The following equipment and reagents are used in the measurements:

- Capel capillary electrophoresis system;
- sodium hydroxide,  $\geq 98$  %;
- sodium tetraborate decahydrate,  $\geq 98$  %;
- sodium cyanoborohydride,  $\geq 95$  %;
- benzocaine (ethyl 4-aminobenzoate),  $\geq 98$  %;
- methanol,  $\geq 99.5$  %;
- chloroform,  $\geq 99$  %;
- acetic acid,  $\geq 99.7$  %;
- D-lactose monohydrate,  $\geq 98$  %.

Data acquisition, collection, processing, and output are performed using a personal computer running under Windows® operating system with Elforun software installed.

### EXAMPLES OF REAL ANALYSES

**BGE:** borate

**Capillary:**  $L_{\text{eff}}/L_{\text{tot}}$  90/100 cm, ID 75  $\mu\text{m}$

**Sample injection:** 250 mbar  $\times$  s


**Voltage:** 25 kV

**Detection:** 305 nm

**Sample:** lactose-free milk

**Found (%):**


1 – lactose ( $< 0.004$ )



**Sample:** lactose-free yoghurt

**Found (%):**

1 – lactose (0.006)



**Sample:** lactose-free yoghurt produced with a violation of technology

**Found (%):**

1 – lactose (0.29)

